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Kirjandus

1. Amstutz U, et al. Clinical Pharmacogenetics Implementation Consortium (CPIC) Guideline for Dihydropyrimidine Dehydrogenase Genotype and Fluoropyrimidine Dosing: 2017 Update. Clinical Pharmacology & Therapeutics 2018;103(2):210-216. [CPIC Guideline for DPYD and Fluoropyrimidines](#)
2. European Medicines Agency. EMA recommendations on DPYD testing prior to treatment with fluorouracil, capecitabine, tegafur and flucytosine. 30. April 2020, EMA/229267/2020. [EMA recommendations on DPD testing prior to treatment with fluorouracil, capecitabine, tegafur and](#)

flucytosine | European Medicines Agency (EMA)

3. LAMP Human DPYD Deficiency KIT (5 variants) LC-DPD5var-LP originaaljuhend, ver 05
4. Pratt VM, et al. DPYD Genotyping Recommendations: A Joint Consensus Recommendation of the Association for Molecular Pathology, American College of Medical Genetics and Genomics, Clinical Pharmacogenetics Implementation Consortium, College of American Pathologists, Dutch Pharmacogenetics Working Group of the Royal Dutch Pharmacists Association, European Society for Pharmacogenomics and Personalized Therapy, Pharmacogenomics Knowledgebase, and Pharmacogene Variation Consortium. The Journal of Molecular Diagnostics, 2024;26(10): 851–863. **DPYD Genotyping Recommendations - The Journal of Molecular Diagnostics**

Koostanud: Krislin Raus, molekulaardiagnostika laborispetsialist

28.09.2026

Viimati uuendatud 29.09.2026